



Netis®



Ansa a rete da recupero

LOOP NET FOR RETRIEVAL

Meditalia®

Netis®



Sezione interna dell'ansa
di forma piatta
*Internal section of
the loop flat shape*



Punta arrotondata atraumatica
Atraumatic rounded tip



Disponibile anche rotante
Available rotatable

Meditalia®

ANSA A RETE DA RECUPERO LOOP NET FOR RETRIEVAL

**Catetere in HDPE con diametro interno maggiorato
essenziale per una migliore scorrevolezza del device**

*HDPE catheter with increased internal diameter
essential for a better scrolling of the device*

**L'ansa monofilamento super stiff aumenta la
stabilità del loop durante il recupero**

*The super stiff monofilament snare increases the
loop stability during retrieval*

**Il filo piatto all'interno che collega il manipolo all'ansa minimizza la frizione durante
la fuoriuscita del catetere**

*Flat wire inside that connects the handle to the snare minimizes friction during
release of the catheter*

Products info

References	Working length	Sheath diameter	Channel	Loop net diameter	Sheath	Morphology
MED-132-NET	230 cm	2,5 mm	2,8 mm	oval 40x65 mm	HDPE	Ultra Stiff
MED-133-NET	230 cm	2,5 mm	2,8 mm	oval 30x60 mm	HDPE	Ultra Stiff
MED-194-NET	160 cm	1,8 mm	2,0 mm	oval 25x45 mm	HDPE	Ultra Stiff
MED-195-NET	230 cm	2,5 mm	2,8 mm	oval 25x45 mm	HDPE	Ultra Stiff
MED-217-NET	260 cm	2,5 mm	2,8 mm	oval 30x60 mm	HDPE	Ultra Stiff
MED-218-NET	230 cm	2,5 mm	2,8 mm	octagonal 40x55 mm	HDPE	Ultra Stiff
MED-219-NET	160 cm	2,5 mm	2,8 mm	octagonal 40x55 mm	HDPE	Ultra Stiff
MED-200-NET	120 cm	1,8 mm	2,0 mm	oval 10x25 mm	HDPE	Ultra Stiff
MED-201-NET	230 cm	2,5 mm	2,8 mm	oval 50x80 mm	HDPE	Ultra Stiff
MED-202-NET	160 cm	2,5 mm	2,8 mm	oval 50x80 mm	HDPE	Ultra Stiff
MED-203-NET	230 cm	2,5 mm	2,8 mm	oval 30x60 mm	HDPE	Ultra Stiff
MED-204-NET	230 cm	2,5 mm	2,8 mm	oval 40x65 mm	HDPE	Ultra Stiff
MED-214-NET	230 cm	2,5 mm	2,8 mm	oval 25x45 mm	HDPE	Flexible
MED-215-NET	230 cm	2,5 mm	2,8 mm	oval 30x60 mm	HDPE	Flexible
MED-216-NET	230 cm	2,5 mm	2,8 mm	oval 40x65 mm	HDPE	Flexible

Meditalia®

www.meditaliasas.it